

PLASMA THAWING BATH

▶ Product Description

LABWAN Plasma Thawing Baths are designed for quick thawing of plasma at 37.0°C with high accuracy. The microprocessor temperature controller maintains accurate temperature conditions with an accuracy of $\pm 0.2^{\circ}\text{C}$ and provides audio-visual alarms in case of temperature deviation. It features a clear-view top cover and a water re-circulation system that ensures uniform temperature. The internal chamber is made of stainless steel 304 for durable and corrosion-free performance. The exterior is made of powder-coated CRCA steel and insulated with a thick layer of Polyurethane Foam (PUF) to maintain temperature and improve energy efficiency. The advanced Plasma Thawing Bath offers a combination of contemporary cabinet design for demanding laboratory environments.

▶ Features

- Microcontroller-based PID controller for superior temperature accuracy.
- Controller displays both SV (Set Value) and PV (Process Value).
- Temperature safety high-low audio-visual alarm.
- High-density, CFC-free Polyurethane Foam (PUF) for minimal temperature loss and reduced power consumption.
- Heavy-gauge stainless steel wire (304) utility basket.
- Water re-circulation system ensures uniform temperature.
- Long-life heating elements for quick response to the controller.
- Attractive front panel bezel for easy cleaning.



Technical	Specification
Product Name	Plasma Thawing Bath
Model name	LW-PTB-01
Operating Voltage	230 VAC
Heating Power	500 W
No. of plasma Bags.	12
Volume	25 Liters
Temp. Range	37.0° C
Temp. accuracy	0.5° C
Stabilization Time	Approximately 30 minutes.
Dimension W x D x H mm	700 W x 370 D x 420 H mm

Application/Uses

- Thawing frozen plasma at 37°C.
- Used in blood banks and hospitals.
- Suitable for clinical and pathology laboratories