



Laboratory Ultrasonic Cleaner



Digital control



Degas



Semiwave



OEM&ODM

Model	LWUC-2	LW-UC3.2	LW-UC4.5	LW-UC6.5	LW-UC10	LW-UC15	LW-UC22	LW-UC30
Capacity	2L/0.53gal	3.2L/0.85gal	4.5L/1.19gal	6.5L/1.72gal	10L/2.64gal	15L/3.96gal	22L/5.81gal	30L/7.93gal
Tank size(mm)	150*135*100	240*135*100	300*150*100	300*150*150	300*240*150	330*300*150	500*300*150	500*300*200
Overall size(mm)	175*160*210	270*165*230	330*175*235	330*175*275	330*270*275	360*325*280	530*330*280	530*330*333
Frequency	40kHz	40kHz	40kHz	40kHz	40kHz	40kHz	40kHz	40kHz
Ultrasonic power	60W	120W/60W	180W/90W	180W/90W	240W/120W	360W/180W	480W/240W	600W/300W
Heating power	100W	100W	200W	200W	300W	300W	500W	500W
Timer	0-30min	0-30min	0-30min	0-30min	0-30min	0-30min	0-30min	0-30min
Temp setting	RT - 80°C	RT - 80°C	RT - 80°C	RT - 80°C	RT - 80°C	RT - 80°C	RT - 80°C	RT - 80°C
Master carton	1PCS	1PCS	1PCS	1PCS	1PCS	1PCS	1PCS	1PCS