

ORBITAL SHAKER

Flat Pan Tray

Product Description

Orbital shaker is a laboratory instrument used to gently mix liquids by moving a platform in a circular (orbital) motion. It is widely used in biology, microbiology, biochemistry, and chemistry laboratories for applications that require uniform mixing without vigorous agitation.

Features

- Maximum 200RPM
- OS20 - LCD Display
- OS21,23,24 - LED Display
- Timer Function available

Accessories

- Main Instrument
- Instruction Manual
- Test Certificate



Technical	Specification
Model	LW-FPT-268
Plate Type	Flat Pan Tray
Plate Size	268 x 268 mm
Speed & Timer	LCD Display
Make	LAWMAN
Max. Load Capacity	2 kg
Speed Range	40 ~ 200 rpm
Orbit Diameter	20 mm
Timer Function	Yes
Ambient Temperature	10°C to 35°C
Power Supply	Adaptor Using AC 220V / 50-60 Hz
Machine Dimension (LxWxH)	295 x 325 x 125 mm
Packaging Dimension (LxWxH)	395 x 370 x 215 mm
Net / Gross Weight	3.1kg / 4.2kg
Packing Type	Cardboard Carton Box
Accessories	Adaptor, Silicon Pad