

ORBITAL SHAKER

Adjustable Clamp Tray

Product Description

Orbital shaker is a laboratory instrument used to gently mix liquids by moving a platform in a circular (orbital) motion. It is widely used in biology, microbiology, biochemistry, and chemistry laboratories for applications that require uniform mixing without vigorous agitation.

Features

- Maximum 200RPM
- OS20 - LCD Display
- OS21,23,24 - LED Display
- Timer Function available

Accessories

- Main Instrument
- Instruction Manual
- Test Certificate



Technical	Specification
Model	LW-OBS-20MM
Plate Type	Adjustable Clamp Tray
Plate Size	250 x 250 mm
Speed & Timer	LCD Display
Make	LAWMAN
Max. Load Capacity	2 kg
Speed Range	40 ~ 200 rpm
Orbit Diameter	20 mm
Timer Function	Yes
Ambient Temperature	10°C to 35°C
Power Supply	AC 220V / 50-60 Hz
Machine Dimension (LxWxH)	295 x 338 x 130 mm
Packaging Dimension (LxWxH)	410 x 300 x 170mm
Net / Gross Weight	4.2kg / 5kg
Packing Type	Cardboard Carton Box
Accessories	Power Cord