

Single Head Vacuum Pressure Pump

Product Description

LABWAN 15 L Single Head Vacuum Pressure Pump is a compact and reliable pump designed for continuous vacuum and pressure operations. It is suitable for laboratory applications such as vacuum filtration, suction, vacuum drying, sample preparation, evaporation and concentration processes, air extraction, and other continuous-duty laboratory operations.

Features

- 1/16 HP
- 220 V
- 1440 RPM
- Vacuum: 22"
- Weight: 3 kg



Technical Specifications

Specification	Details
Product Name	15 L Single Head Vacuum Pressure Pump
Horse Power	1/16 HP
Power Supply	220 V
Motor Speed	1440 RPM
Vacuum	22"
Weight	3 kg
Maximum Vacuum	13 mbar
Work Mode	Continuous Operation
Nozzle Size	8 mm
Power	95 W
Brand	LABWAN