

LABWAN ROLLER MIXER

Product Description

The LABWAN LWRM-6 Roller Mixer provides gentle and uniform mixing through rocking and rolling motion, with digital speed and timer control. It features six rollers of 28 cm length, a speed range of 20–80 RPM, a timing range of 1 minute to 23 hours 59 minutes, and compatibility with tube diameters of 6–34 mm. It supports up to 10 × 50 ml centrifuge tubes or 0.5 kg load and operates on AC 220 V / 50 Hz with 12 W power consumption.

Features

- Rocking and rolling motion for complete mixing
- Digital control of time and speed
- Six convenient rollers
- Roller length: 28 cm
- Suitable for various tube diameters
- Stable functioning and very low noise

Accessories

- Adapter
- Instruction Manual



Parameter	Specification
Product name	LABWAN Roller Mixer
Model name	LWRM-6
Make	LABWAN
Motion	Rocking and Rolling
Speed Range	20 to 80 RPM
Rollers	6 rollers
Max Capacity	10 * 50ml Centrifuge tube / 0.5kg
Timing Range	1min to 23h 59mins
Compatible Tube Dia	Ø6 to 34mm
Roller Length	28 cm
Power	12W
Power Supply	AC 220V / 50 Hz
Dimensions & Weight	445 × 260 × 103 mm & 3 kg